

## Lipids External Quality Control Program

### Directions for performing test:

- Select the appropriate sample as indicated by the label on the vial
- Carefully add exactly 0.5 mL of distilled water (or water for injection) and allow the contents to dissolve for 15 minutes.
- Swirl the contents until homogenous
- Run the sample as a Patient test with the patient ID being the sample number indicated in the bottle
- Draw the solution into the plastic pipette and apply a drop of control solution on to the disc's suction point, ensuring the sample area is completely filled.
- Press the hinge cover down firmly until two clicks are heard and place the disc into the machine.
- Once analysis is complete, fill out this sheet and fax results to the iCCnet SA office (08) 7117 0635

SAMPLE NUMBER

Total Cholesterol (mmol/L)

LDL Result (mmol/L)

HDL Result (mmol/L)

Triglycerides Result (mmol/L)

Date/Time Reconstituted

Date/Time Analysed

Operator



Site Name

**Please fax results to iCCnet SA on 08 7117 0635**