

StatStrip Glucose/Ketone

External Quality Control Program

Directions for performing test:

- Select the appropriate sample as indicated on the pipette (blue solution) and cut off the tip.
- Run the control solution as a patient test, scanning the sample barcode as the patient ID.
- Discard the first few drops and apply a large drop to the test strip ensuring there are no air bubbles.
- Test all StatStrip meters using the same solution.
- Once analysis is complete, fill out this form and fax to iCCnet on (08) 7117 0635 or email Health.iCCnetSupport@sa.gov.au.

		EQA:		EQA:	
	Device Serial Number	Glucose mmol/L	Ketone mmol/L	Glucose mmol/L	Ketone mmol/L
Device 1					
Device 2					
Device 3					
Device 4					

Date Analysed:

Operator:

Site Name:

Please fax results to iCCnet SA on 08 7117 0635